IKA

designed for scientists

IKA RV 10 control

IKA RV 10 auto

Operating instructions	EN	4		
Veiligheidsaanwijzingen	NL	28	Bezpečnostné pokyny	SK 44
Säkerhetsanvisningar	SV	30	Ohutusjuhised	ET 46
Sikkerhedshenvisninger	DA	32	Drošības norādījumi	LV 48
Sikkerhedsveiledning	NO	34	Saugos nurodymai	LT 50
Turvallisuusohjeet	FI	36	Указания за безопасност	BG 52
Bezpečnostní pokyny	CS	38	Instrucțiuni de siguranță	RO 54
Biztonsági utasítások	HU	40	Υποδείξεις ασφαλείας	EL 56
Varnostna navodila	SL	42		

- EN**
- A RV 10 auto drive
 - B Heating bath
 - C Bracket
 - D Woulff bottle
 - E Glassware

- EN**
- A "Power" switch
 - B "Power" key
 - C Lift position "▲" key
 - D Lift position "▼" key
 - E "Rotating speed" key
 - F "Vent" key
 - G "Vacuum" key
 - H "Menu" key
 - J "Back" key
 - K Start/Stop knob

Contents

	Page		
Device setup	2	Interfaces and outputs	23
EU Declaration of conformity	3	Maintenance and cleaning	24
Explication of warning symbols	3	Accessories	24
Safety instructions	4	Error messages	24
Intended use	5	Warranty	26
Unpacking	6	Technical data	27
Useful information	6	Drawings for mounting the glassware	58
Vacuum control	7	Drawings for connecting the hose system	60
Setting up	8	Solvent table (excerpt)	63
Commissioning	14		

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at sales@ika.com.

Explication of warning symbols

Indicates an imminently hazardous situation, which, if not avoided, will result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

Indicates practices which, if not avoided, can result in equipment damage.

Danger - Reference to exposure to a hot surface!

Safety instructions

For your protection

- Read the operating instructions completely before starting up and follow the safety instructions.
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations. When working under a vacuum in particular!
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There may be a risk of:
 - splashing liquids,
 - body parts, hair, clothing and jewellery getting caught,
 - injury as a result of glass breakage.

DANGER

Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapors, dusts or biological and microbiological materials can be hazardous to user.

- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
 - stress due to incorrect assembly,
 - external mechanical hazards,
 - local temperature peaks.
- Ensure that the stand does not start to move due to vibrations respectively unbalance.
- Beware of hazards due to:
 - flammable materials,
 - combustible media with a low boiling temperature.

CAUTION

Only process and heat media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set. The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used. (acc. to EN 61010-2-010)

- Do not operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Tasks with the device must only be performed when operation is monitored.
- Operation with excess pressure is not permitted (for cooling water pressure see "Technical Data").
- Do not cover the ventilation slots of the device in order to ensure adequate cooling of the drive.
- There may be electrostatic discharges between the medium and the drive which could pose a direct danger.
- The appliance is not suitable for manual operation (except lift movement).
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.

- Refer to the operating instructions for the accessories, e.g. vacuum pump, heating bath.
- Position the positive pressure outlet of the vacuum pump under a fume hood.
- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Adapt the quantity and the type of distill and to the size of the distillation equipment. The condenser must work properly. Monitor the cooling water flow rate at the condenser outlet.
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at condenser) in order to prevent a pressure build-up.
- Please note that dangerous concentrations of gases, vapors or particulate matter can escape through the outlet at the condenser. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.
- Evacuated glass vessels must not be heated only on one side; the evaporating flask must rotate during the heating phase.
- The glassware is designed for operation under a vacuum of up to 1 mbar. The equipment must be evacuated prior to heating (see chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, uncondensed vapors must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.

CAUTION

Avoid peroxide formation. Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing! Keep liquids that tend to form organic peroxides away from light, in particular from UV rays and check them prior to distillation and exhaust for the presence of peroxides. Any existing peroxides must be eliminated. Many organic compounds are prone to the formation of peroxides, e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.

DANGER

The heating bath, tempering medium, evaporation flask and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.

WARNING

Avoid delayed boiling! Never heat the evaporating flask in the heating bath without switching on the rotary drive! Sudden foaming or exhaust gases indicate that the flask content is beginning to decompose. Switch off heating immediately. Use the lifting mechanism to lift the evaporation flask out of the heating bath. Evacuate the danger zone and warn those in the surrounding area!

Safety lift

When the device is switched off or the power supply disconnected, the internal safety lift removes the evaporating flask from the heating bath.

The safety lift at loss of power is designed for a maximum total weight (glassware and solvent) of 3.1 kg.

Sample calculation of the maximum load with vertical glassware and a 1 litre flask:

Condenser + receiving flask + evaporating flask + fittings =

1200 gr + 400 gr + 280 gr + 100 gr = 1980 gr

Maximum loading of solvent = 3100 gr – 1980 gr = 1120 gr

Due to the design, operation of the safety lift cannot be guaranteed for higher loads!

When using other types of condensers such as dry ice or intensive condensers, and also when using return distillation distributors with slip-on condensers, it may be necessary to reduce the load by the amount of the added weight of the glass apparatus.

Thus, prior to distillation, check whether the lift goes up without power when laden with the glassware and distillation material.

The safety lift must be checked daily prior to operation. Use the motor to move the lift to the lowest position and press the "Power" key on the front plate or the main switch on the back right side of the device using the maximum total weight of 3.1 kg.
→ The evaporating flask is lifted out of the heating bath.

When using for the first time or after a long period of inactivity, proceed as follows: use the motor to move the lift several times to the top and bottom positions by pressing the lift function keys prior to interrupting the power supply.

If the safety lift does not work, please contact the **IKA** Service department.

For the evaporating equipment (evaporating flask plus contents), the maximum permissible weight is 3.0 kg! Loads greater than this risk breakage of glass at the vapor tube!

Ensure that the safety lift has been powered down.

When working with large loads, always use low speeds. Unbalanced loads can result in breakage of the vapor tube!

- A vacuum may be formed inside the glassware in the case of power outage. The glassware must be vented manually.

CAUTION

Never operate the device when the evaporation flask is rotating and the lift is raised. Always lower the evaporation flask into the heating bath first before starting the rotation drive. Otherwise hot tempering medium may be sprayed out!

- Set the speed of the drive so no tempering medium is sprayed out as a result of the evaporation flask rotating in the heating bath. If necessary reduce the speed.
- Do not touch rotating parts during operation.
- Imbalance may result in uncontrolled resonance behavior of the device or assembly. Glass apparatus may be damaged or destroyed. In the event of unbalance or unusual noises, switch off the appliance immediately or reduce the speed.
- The appliance does not start up again automatically following a cut in the power supply.
- The device can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Always check evaporation processes before starting automatic operation. Automatic operation must not be used for unfamiliar evaporation processes. The following automatic modes are available in addition to manual mode:
 - 100% distillation,
 - Volume-/quantity-dependent distillation.

For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids, etc.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts.

Intended use

• Use

Together with the accessories recommended by **IKA**, the device is suitable for:

- quick and gentle distillation of liquids,
- evaporation of solutions and suspensions,
- crystallization, synthesis or cleaning of fine chemicals,
- drying of powder and granulate material,
- recycling of solvents.

Mode of operation: Tabletop device

• Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

• Further remarks

Legal bases and country-specific characteristics, especially for medicine and food:

The operator is responsible for compliance with the laws and regulations applicable to him.


Unpacking

• Unpacking

- Unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post, rail or logistics company).

• Scope of delivery

	Drive RV 10 auto	Heating bath HB digital	Vertical glassware RV 10.1	Vertical glassware RV 10.10 coated	Condensate flask 1	Bracket 2	Cooler locking device, compl. 3	Vacuum hose (2 x 0.55 m) 4	Water discharge hose (1 x 1 m) 5	Vapour tube 6	Ring spanner 7	Handhold 8	RV 10.4002 vacuum valve for laboratory vacuum system 9	IKA Vacstar digital vacuum pump 10	Operating instructions
RV 10 control V	x	x	x		x	x	x	x	x	x	x	x	x		x
RV 10 control V-C	x	x		x	x	x	x	x	x	x	x	x	x		x
RV 10 control FLEX	x	x			x	x	x	x	x	x	x	x	x		x
RV 10 auto V	x	x	x		x	x	x	x	x	x	x	x			x
RV 10 auto V-C	x	x		x	x	x	x	x	x	x	x	x			x
RV 10 auto FLEX	x	x			x	x	x	x	x	x	x	x			x
RV 10 auto pro V	x	x	x		x	x	x	x	x	x	x	x		x	x
RV 10 auto pro V-C	x	x		x	x	x	x	x	x	x	x	x		x	x
RV 10 auto pro FLEX	x	x			x	x	x	x	x	x	x	x		x	x

Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.

The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C). Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation cooler should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a condensate flask and/or a vacuum separator. Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems.

Speed, temperature, piston size and system pressure all affect the evaporator air capacity.

The optimum capacity of the flow-through condenser is approx. 60%.

This corresponds to condensation on approx. 2/3 of the cooling coil.

With larger capacities there is the risk that the uncondensed solvent vapor will be extracted.

In order to prevent this, the condenser loading safety monitoring system can be activated. Before activating the safety monitoring, first carry out several distillations, noting the value given under menu option "EVAPORATOR → Advanced → Max power of cooler". Now calculate the heat rejection rate using the following formula:

$$\text{Cooler power limit} = \text{Max power of cooler} * 1.3$$

Enter this value under menu option "EVAPORATOR → Advanced → Cooler power limit". The device will now stop the measurement and issue an error message if the threshold for the condenser is exceeded. If the configuration of the condenser equipment or the distillation process is modified, then it may be necessary to re-evaluate the maximum condenser loading.

The device is equipped with a piston-operated safety mechanism. The glass apparatus may contain a vacuum following interruption to the power supply; vent the system before switching back on. If the power cuts out, the evaporation flask is automatically lifted out of the heating bath by an integrated gas spring.

CAUTION

The safety lift must be checked daily prior to operation. For more information see the section "Safety Instructions – Safety Lift"!

In addition to offering a full range of manual and semi-automatic evaporation operating modes, the RV 10 control rotary evaporator can also be used for fully automatic and quantity-controlled evaporation processes. Standard fittings for fully-automatic operation include a vacuum controller, a cooling water differential temperature measurement sensor and a cooling water flow rate meter. The vacuum can be controlled in both two-position and speed-controlled pumping operation with the help of the integrated vacuum controller. The device is designed for operation with a cooling water supply system (e.g. laboratory thermostat), but can also be run off a water supply line. Please refer to the Technical Data for information on cooling water pressure, temperature constancy, flow rate and local regulations.

We recommend using the RV 10.5001 water regulator valve when operating directly off a water supply line.

This valve allows the cooling water flow rate to be adjusted and automatically cuts off the cooling water supply on completion of the distillation process.

The RV 10.5002 water filter can be used to prevent particles of dirt passing from the water supply line into the valve mechanism. If pressure reduction is required, we recommend fitting the RV 10.5003 pressure regulator valve in the pipe immediately downstream of the pick-off point.

Automatic volume-dependent distillation: The apparatus must be brought up to operating temperature. This is achieved by conducting a test distillation.

RV 10 auto: The rotary evaporator is fitted as standard with speed control pump, internal venting valve and a pressure measuring chamber.

RV 10 control (for vacuum): The external two-position vacuum control is also available for the application which in-house vacuum system is available. Vacuum valve RV 10.4002 must be used.

Vacuum control

RV 10 control

A vacuum is made in the glassware with the help of a vacuum pump. The vacuum pump operates at constant speed which generally cannot be adjusted.

Once the target value is reached the suction line is interrupted. The device will fall slightly short of the target value setting due to the minimum time difference from detection of the pressure value, the target value comparison for switching the vacuum valve and the suction line of the vacuum generating system. If the pressure in the system increases again due to the natural leakage rate, the valve opens the suction line.


The switching frequency and accuracy can be set via the hysteresis setting.

With the VC 2.4 pump control connected, the power supply to the pump is also switched on and off accordingly. This reduces the noise level and increases the service life of the pump membranes and motor.

The actual value fluctuates in size of the set hysteresis around the target value.

Automatic boiling point recognition is not accurately possible when using two-position control.

Schematic view of a two-position vacuum control

RV 10 auto

Imprecision of the two-position control is avoided by using the speed-vacuum control.

However a vacuum pump with a controlled speed is needed for speed-vacuum control; the RV 10.4002 vacuum valve is not required and must not be connected! To do so, connect the vacuum pump (e.g. IKA Vacstar digital pump) to the port on the rear of the rotary evaporator.

With this kind of control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value.

Once the target value is reached, the pump only operates according to the leakage rate.


It is possible to achieve quieter operation and more exact vacuum control.

The speed-vacuum control on the RV 10 is controlled automatically as soon as an appropriate vacuum pump is attached.

Automatic boiling point recognition is possible with this type of control, i.e. the system reaches and keeps the solvent at its boiling point in automatic mode.

As it is not necessary to know the boiling point of the solvent when in automatic boiling point recognition mode (unlike in the case of volume-dependent distillation), these two operating modes are mutually exclusive.

Schematic view of speed-vacuum control


Setting up

Drive RV 10 auto

Loosen transportation lock!


 **CAUTION**

- Hold the lift with your hand on the height position and remove the thumb screw (1) by turning counterclockwise on the back of the appliance.
- Once the transportation lock has been removed, the lift moves slowly to its upper end position. The distance is approx. 140 mm.
- Connect the supplied power cable to the connection socket (2).

Adjustable base

Note: Only one base foot is adjustable.

Fix the bracket


Insert the bottle and attach the supplied hose connectors to the bottle

Set the drive at an angle of approx. 30°

- Remove the clamping device for the angle setting of the rotation drive on the right side of the lift by rotating the knurled screw counterclockwise (gently press and rotate the knurled screw (1) at the same time to move it further in or out).
- Set the drive at an angle of approx. 30°.

- Then secure the rotation drive from being accidentally turned by tightening the knurled screw in a clockwise direction.

Heating bath

CAUTION

Refer to the chapter "Operation" in the heating bath instruction manual!

- Place the heating bath on the stand of the rotation drive and push it into the left position.



Note: Data is exchanged between the drive unit and the heating bath by means of an infrared link (1). Please note that reliable communication is only guaranteed when the infrared beam has clear line-of-sight to the detector.


Glassware

NOTE: Please read the operating instructions of the glassware for the safe handling of laboratory glassware !

- There is a locking knob for locking or unlocking the locking device.
If a red mark is shown, the locking device is unlocked, otherwise it is locked.
To lock or unlock the locking device, push the locking knob to the end position.

- Open the locking device on the drive head by turning it 60° counterclockwise. Then you can see a red mark .
- Feed the vapour tube in until it stops.
- Then, lock the locking device by turning it clockwise by 60°.
- Pushing the lock knob to the end position, the red mark should be covered and invisible.
- The vapour tube is not allowed to be pulled out!
- Check the correct axial locking device on the vapour tube.
- Keep the red mark invisible.

Fit the condenser seal




- Place the condenser seal in the condenser receptacle and fit the glassware to the device according to the assembly instructions.

First use

- Insert vapor tube (1).
- Insert the vacuum seal (2). Pay attention to the correct position of the vacuum seal!
- Slide the union nut (3b) over the flange on the condenser (3a).
- Also slide the annular spring (3c) over the flange on the condenser (3a).
- Position the condenser (3a) on top of the seal (2).
- Tighten the union nut (3b) by hand onto the threaded flange. Run the device for 20 minutes at 120 rpm. Finally, retighten the union nut (3b) by hand.

Note: Follow the mounting instructions for the glassware.

No stand is required if the condenser has been properly fitted and secured by tightening the condenser union nut on the rotary drive.

The purpose of the stand is solely to prevent the condenser from twisting.

CAUTION

If the stand is fitted incorrectly, the strong forces at the Velcro strap can cause stresses within the glass, which in turn can cause damage to the glass condenser.

After the condenser has been fitted, the stand is attached to the rotary drive. Make sure that the condenser is installed parallel to the lift body.

Assembling the vertical glassware cooler locking device

- Assemble the condenser locking device according to the diagram.
- Mount the plate (1) using the thumb screw (2).
- Put the support rod (4) on the plate (1) and attach it with the screw nut (3).
- Attach the rubber protector (5).
- Fasten the Velcro® band (6) to the support rod (4).
- Secure the vertical glassware with the Velcro® band (6).

Mounting the glassware

Note: Observe "Drawings for mounting the glassware" on page 58.

Item	Designation	Quantity					
		RV 10.1 non-coated RV 10.10 coated	RV 10.2 non-coated RV 10.20 coated	RV 10.3 non-coated RV 10.30 coated	RV 10.4 non-coated RV 10.40 coated	RV 10.5 non-coated RV 10.50 coated	RV 10.6 non-coated RV 10.60 coated
1	Receiving flask	1	1	1	1	1	1
2	Clamp (stainless steel)	1	1	1	1	1	1
4	Stopcock	1	1	1	1	1	1
5	Pipe tube	1	-	1	1	1	1
6	Condenser	1 Vertical condenser	1 Diagonal condenser	1 Vertical-intensive condenser	1 Dry ice condenser	1 Vertical condenser	1 Vertical-intensive condenser
7	Connection	1 Vacuum connection	1 Introduction sleeve	1 Vacuum connection	-	1 Vacuum connection	1 Vacuum connection
8	Clamp (plastic)	1	1	1	-	1	1
10	Evaporation flask 1000 ml	1	1	1	1	1	1
11	Ball joint clamp RV 05.10	1	1	1	1	1	1
12	Screw joint cap	4	4	4	2	4	4
13	Hose connection	4	4	4	2	4	4
14	Manifold	-	-	1	-	1	1
15	Cover cap	-	-	-	1	-	-
16	Cartridge	-	-	-	1	-	-
17	Ring, slotted	-	-	-	1	-	-
18	Washer	1	-	1	1	1	1
19	O-ring	-	-	-	1	-	-
20	Plug, high vacuum	-	-	-	-	1	1


Note:

The PTFE tube (5) with the washer (18) can be mounted optionally on the vertical glass sets.


It serves for filling the evaporator flask when there is a vacuum in the glass set.

Solvent can be sucked through the PTFE tube into the evaporator flask by opening the stopcock (4).

Mounting the washer

Drainage washer

Note: Pay attention to the correct position of the washer.


Removing the condenser

- Use the ring spanner provided to loosen union nuts that are tightly fitted.
- Loosen the union nut by turning counterclockwise.
- Remove the Velcro®.

Note: Ring spanner is only used for removing the condenser, while fixing the condenser please use hands.

Hose system

Note: Observe "Drawings for connecting the hose system" on page 60.


RV 10 control operation

- Fit RV 10.4002 vacuum, in the bracket provided for this purpose and connect the vacuum hose to the valve.


RV 10 auto operation

- No additional valve is required for vacuum normal operation with a speed-vacuum control pump. Connect the pump's vacuum inlet directly to the RV 10 pump hose connection.
- Connect the pump controller cable to RV 10 control.

Connect RV 10.4002


Note: When connecting the pump controller cable, make sure that the snap-and-lock plug has been inserted into the stop position inside the pump controller cable connection.

Connect RV 10 temperature sensor (dT)


- Connect the temperature sensor to the lower socket (dT). Always check for correct direction of the arrow.

Note: When connecting the temperature sensor, make sure that the snap-and-lock plug has been inserted into the stop position inside the temperature sensor connection.


Water

- Connect the water inlet to the water supply (with optional RV 10.5001 water regulator valve). Please refer to the technical data relating to the water supply. The RV 10.5001 water regulator valve is not suitable for use on the condenser unit because it would reduce the flow rate too greatly.
- Connect the water hoses to the glass condenser (short hose (1) = lower outlet, long hose (2) = upper inlet) and fix hoses in place.


- Connect the water outlet hose (supplied) by pushing it onto the nipple until the stop inside the connection socket is reached.
Note: Ensure that the condenser inlet and condenser outlet are connected correctly. Automatic distillation cannot be carried out if the water inlet and outlet hoses are not fitted according to the instructions and illustrations, since this would cause the temperature readings for the inlet and outlet to be measured incorrectly.


- Release the hose connection using the supplied tool.

- Insert the valve connector cable (RV 10.5001, RV 10.4002) into the appropriate socket. Connect the speed-vacuum control pump. As a result, the RV 10 control automatically switches to normal speed-vacuum control operation.

- Connect "LOAD" to the Condensate flask and the Condensate flask to the vacuum connector on the condenser using the vacuum hoses supplied. Please note that the vacuum must always be connected at the highest connection point on the condenser.


Commissioning

Working screen at the time of delivery

The start screen appears for a few seconds after the device is switched on. The device name and the software versions are displayed.

After this, the working screen appears automatically on the display.

An information prompt then appears to download the **IKA** firmware update tool.

Explanation of symbols on the working screen

The symbols displayed change depending on the status and settings of RV 10 auto.

Symbol	Designation	Description
	Bluetooth®	This symbol means the RV 10 is communicating via Bluetooth® with a PC or tablet. The symbol no longer appears if no Bluetooth® communication is being performed.
	USB	This symbol means RV 10 is communicating via a USB cable. The symbol no longer appears if no USB cable is being used for communicating with the station.
LIMIT	Limit	This symbol indicates if a minimum or maximum limit of a set value is reached.
	Temperature sensor	This symbol appears when the temperature display is active in the display.
PC	PC control	This symbol means that device is connected to a computer and that it is being controlled by the computer.
PR	Program controlled	This symbol means that device is being controlled by a program.
	Rotation	This symbol indicates that the rotation is in continuous mode and indicates the direction of rotation.
	Drying	This symbol indicates that the drying option was selected (only in 100 % mode).
	Flow rate	This symbol indicates that a fluid (water) is running through the system.
	Cooler connected	This symbol indicates that a cooler is connected and running.
	Venting	This symbol indicates that the venting valve is open. If this symbol is not shown the venting valve is closed.
vac	Vacuum	This symbol indicates that the pressure inside the RV 10 system is lower than the atmosphere pressure.
	Detecting	This symbol is shown during the automatic evacuation mode.
	Boiling	This symbol indicates that a boiling point was found (only in automatic mode).

Menu structure

			Factory settings	
Evaporator	Modes	Automatic	Bath medium.....	Water
			Start temperature.....	60 °C
		Manual.....		1013 mbar
		Pump %.....		50%
		Volume	Solvent.....	Acetic acid
			Target.....	100 ml
			Efficiency.....	80%
			Adjustment.....	-
		Program.....		-
		100%	Solvent.....	Acetic acid
	Drying.....		-	
	Cleaning	Start after...	mm:ss.....	30:00
		Duration	mm:ss.....	06:00
		Pump speed	Value.....	20%
		Clean now.....		-
	Limits	Manual	Minimum.....	2 mbar
			Maximum.....	1013 mbar
		Pump %	Minimum.....	0%
			Maximum.....	100%
	Hysteresis	Automatic.....		✓
		Manual	Value.....	±15 mbar
	Solvent library	Add solvent		
		Add solvent		
		Add solvent	Name.....	-
		Add solvent	Formula.....	-
		Add solvent	Bath temperature.....	-
			Rotation speed.....	-
			Boiling point.....	-
			Efficiency.....	-
			Heat capacity.....	-
			Enthalpy.....	-
			Density.....	-
			i-factor.....	-
			Acetic acid	
			Acetone	
		Acetonitrile		
	...			
	Xylene			
	Name.....	-		
	Formula.....	-		
	Bath temperature.....	-		
	Rotation speed.....	-		
	Boiling point.....	-		
	Efficiency.....	-		
	Heat capacity.....	-		
	Enthalpy.....	-		
	Density.....	-		
	i-factor.....	-		
Calculator	Name.....		Acetic acid	
	Vapor temperature		40 °C	
	Heating bath temperature.....		60 °C	
	Pressure.....		46 mbar	
Actions after distillation	Stop heating.....		-	
	Stop rotation.....		✓	
	Lift up.....		✓	
	Close water valve.....		✓	
	Open venting valve.....		✓	
	Cleaning.....		-	

Menu


		Factory settings				
	Advanced	Flow control.....	✓			
		Bath medium.....	Water			
		Cooler power limit.....	900 w			
		Max power of cooler.....	0 w			
	Rotation	Interval mode	Interval	Ccw Value.....	00:10	
				Cw Value.....	00:10	
		Speed limit	Activate.....	-		
			Minimum.....	5 rpm		
			Maximum.....	300 rpm		
	Timer	Timer function	Time	hh:mm:ss.....	00:00:00	
			Beep after timeout.....	-		
			Activate.....	-		
		Display.....	-			
	Display	Timer.....	✓			
		Flow rate.....	✓			
		Δ T.....	✓			
		Cooler power.....	✓			
	Programs	Program 1-10	Select.....	-		
			Edit	Edit	Pressure.....	-
					Rotation (Speed).....	-
					hh:mm:ss.....	-
			Insert.....	-		
		Delete.....	-			
		Last Measurement	Save as.....	-		
			Edit	Edit	Pressure.....	-
					Rotation (Speed).....	-
					hh:mm:ss.....	-
Insert.....	-					
Delete.....	-					
	Safety	Password	Value.....	0		
			Value.....	0		
			Value.....	0		
		Continue after power failure	Activate.....	-		
			Time	mm:ss.....	00:10	

		Factory settings		
 Service	Valves	Venting valve.....	-	
		Vacuum valve.....	-	
		Water valve.....	-	
	Pump.....	-	-	
	Temperature adjustment.....	-	-	
 Settings	Language	English.....	✓	
		Deutsch.....	-	
		Français.....	-	
		Español.....	-	
		Italiano.....	-	
		Português.....	-	
		Русский язык.....	-	
		中文.....	-	
		한국어.....	-	
		日本語.....	-	
	Units	Pressure	mbar.....	✓
			hPa.....	-
			mmHg.....	-
			Torr.....	-
		Amount	mL.....	✓
			g.....	-
	Display	Background	Black.....	✓
			White.....	-
		Brightness Value.....	100%	
		Firmware update info.....	-	
	Sound	Volume Value.....	100%	
		Key Tone.....	-	
	Factory settings	Menu values.....	-	
		Programs.....	-	
		Solvent library.....	-	
		All.....	-	
	Communication	Device name.....	RV 10 auto	
		Bluetooth.....	✓	
		Labworldsoft 5 protocol.....	✓	
	Information	Firmware update info.....	-	
		Display version.....	-	
		Logic version.....	-	
Pressure max.....		1100 mbar		
Pressure min.....		1 mbar		
Pump % max.....		100%		
Pump % min.....		0%		
Rot. speed max.....		300 rpm		
Rot. speed min.....		5 rpm		
Operating hours.....		0 hours		

Change the menu items, values or settings by turning the rotary knob and confirm by pushing it.

1. Evaporator

Modes

Automatic: In this mode, the system will detect the actual evaporation by checking the temperature difference between water input and water output. The default setting of the medium temperature is 60 °C, the user can set the medium temperature manually. If the expected evaporation temperature is higher than 90 °C, please change the medium from water to oil, and use high temperature medium, e.g. silicone oil. In this mode, measurement will start automatically when the heating bath reaches the set temperature, the temperature difference will be below a threshold value, the water flow rate must be within 30-100 L/h, and the lift moves down automatically. You can change the speed manually.

Manual: In this mode the user can set all values manually. When the measurement is started the system is evacuated until it reaches the set pressure value. The lift has to be moved manually.

Pump %: In this mode the pump can be run in a long term, setting a value between 100 % and 1% of the possible pump speed.

Volume: This mode is used to distill a specific amount of the used solvent. The measured values of the cooling water flow rate and the cooling water temperature difference from the basis for a heat balance calculation used to determine the quantity of distillate at every stage of the distillation. The distillation process is stopped when the specified quantity of distillate has been reached.

Solvent: Select the solvent you want to distill.

Target: Enter the amount, which you want to distill off.

Efficiency: Balancing requires a precise thermal efficiency to be entered.

Adjustment: The efficiency is dependent on a combination of factors that affect the distillation process, so it is initially estimated. The first time a distillation process is carried out should therefore be used for adjustment. The actual efficiency will be recalculated after the measurement based on the target volume and the actually distilled volume.


Program: This entry is a link to the "Programs" menu.

100 %:

Solvent: The distillation process is stopped when the measured cooling water temperature difference falls below a threshold value, i.e. as soon as a solvent has been completely distilled off.

Drying: The cooling water temperature difference is not monitored if this option is activated, e.g. use in processes for drying powdered media.

Cleaning

Start after...: Set the minimum time that a measurement has to run before the cleaning after the measurement is activated. The cleaning process starts when the measurement has been running for at least the set time.

Duration: Set the duration of the cleaning process.


Pump speed: The set pump speed in the cleaning menu is used to set the speed of the pump during the cleaning process.

Higher speed causes partial vacuum in the glassware if the system is closed. High speed cleaning affects the pump membrane lifetime.

Clean now: This starts the cleaning process manually.


Limits

Manual: This option allows user to set the permitting operating vacuum range.

Pump %: This option allows user to set the permitting operating pump speed range.

Hysteresis


The (vacuum) hysteresis value describes the difference between the closing and opening pressures for the vacuum valve. Setting the hysteresis value too low can have a negative effect on the vacuum pump and the vacuum valve. The hysteresis function is only used for indoor vacuum or uncontrollable vacuum pumps.

Automatic: Hysteresis is always 10 % of the actual pressure.

Manual: Manual value specification of hysteresis.

Solvent Library


User can select the solvent which will be distilled, the parameters (such as, Name, Formula, Bath temperature, Rotation speed, Boiling point, Efficiency, Heat capacity, Enthalpy, Density and i-factor) will be shown in the library. User can set the evaporating bath temperature, efficiency and the rotation speed.

It is possible to add up to five new customized solvents to the solvent list. The customized solvents appear in orange and can be found on the top of the solvent library.


Name and formula: The name and the formula help to identify the solvent.

Calculator


It helps you to calculate the boiling point (vapor temperature) of solvent, heating bath temperature and pressure. If one of the 3 values has been defined, the other 2 values will be given.

Actions after distillation

It is possible to define actions when the measurement is stopped (manually or automatically).

Advanced


Flow control: This setting can be deactivated in the manual mode if dry ice is used instead of a cooler.

Bath medium: Select either water or oil as the medium for the heating bath. If water is the medium, please select water, the temperature will be from 20 °C to 90 °C. If oil is the medium, please select oil, the temperature range will be from 20 °C to 180 °C.

Cooler power limit: This setting is used to limit the cooler power. If the actual cooler power is higher than the cooler power limit a warning appears and the measurement is stopped.

Max power of cooler: This value shows the maximum cooler power value of a measurement. This value serves as information and cannot be changed.

2. Rotation

Interval mode

The rotation can be run in an interval mode of alternating the rotation direction.

Interval: The interval sets the time how long the rotation takes place in one direction.

Activate: This entry activates the interval mode for the rotation.

Speed limit

The "Speed limit" menu allows the user to set the desired maximum and minimum speed limit for the rotation.

3. Timer

Timer function

Time: This setting allows the user to specify the actual time for the measurement procedure in the manual mode. A default time can also be set for the timer, which allows the user to start the measurement task for a standard time. The device stops automatically after expiry of the set time, and the set time used for the measurement procedure appears on the display.

Note: The user can stop the measurement procedure before expiry of the set time. In this case the countdown of the timer is interrupted.

Beep after timeout: This menu item allows the user to turn on/off the sound that occurs when the timer reaches 00:00:00. A tick shows that the option is activated.

Activate: In this menu the user can activate or deactivate the timer function. A tick shows that the option is activated.

Display


In the "Timer" menu the user can specify that the timer is displayed on the display/working screen. A tick shows that the option is activated.

4. Display

In this menu the user can set which values (Timer, Flow rate, ΔT and/or Cooler power) should be shown on the main screen. A tick shows that the option is activated.

5. Program

Program 1 - 10

Under menu "Programs", 10 user-defined pressure-rotation (rpm)-time profiles can be created. In addition, the user can define whether the interval mode is activated or not in the programs. If the interval mode is activated, the run time/stop time value from "Interval Mode" setting will be taken.

Note: If the user needs to activate the interval mode in one segment of a program, he should set the "Ccw/cw time" in menu option "Interval", meanwhile activate the "Ccw/cw" function in menu option "Interval Mode" (see section "Rotation").

Select: Select the program.

Edit: Edit the selected program parameters. Start to edit the selected program parameters by pressing menu option "Edit" with start/stop knob. The user can edit, insert or delete one selected program segment in the program. When user edits the program time for at least one segment, a tick (✓) for respective program will appear.

Delete: Delete the selected program. If a selected program is deleted by pressing on menu option "Delete" with start/stop knob, all the program parameters will be emptied. The tick (✓) disappears.

Last measurement

Save as: Save the distillation sequence as a program.

Edit: Edit the selected program parameters.

Start to edit the selected program parameters by pressing on menu option "Edit" with start/stop knob. The user can edit, insert or delete one selected program segment in the program. When the user edits the program time for at least one segment, a tick (✓) for the respective program will appear.

Delete: Delete the selected program. If a selected program is deleted by pressing on menu option "Delete" with start/stop knob, all the program parameters will be emptied. The tick (✓) disappears.

Details for editing the program

When a program is being edited, the following screen appears.

No.	Pressure	Rotation	h:m:ss
01	500 mbar	40 rpm	00:00:20
02	400 mbar	50 rpm	00:00:20
03	300 mbar	60 rpm	00:00:20
04	200 mbar	70 rpm	00:00:20
05	100 mbar	80 rpm	00:00:20

edit insert delete


In this program, user can define up to 10 segments. The selected segment is highlighted. Then, the user can edit, insert or delete a segment in this program. When the "Back" button is pressed after editing, the program is saved automatically.

Edit: When the background of the selected value turns yellow the user can change the pressure, the speed or the time.


Insert: With the insert option, a new segment will be insert below the selected segment.

Delete: When deleting a highlighted segment, the settings will be emptied. The section jumps to next segment.


Sample for editing the program

Example of saving the last measurement

6. Safety

Password

It is possible to set a password to secure the main menu from setting changes. To deactivate the password set the password 0 0 0.


Continue after power failure

If this option is activated the measurement will continue after a break caused by a power failure. This option is only available in Automatic, 100 % and Volume mode.

Activate: If tick is shown the option is activated.

Time: The measurement will be continued if the the power comes back within the set time.

7. Service

To check functionality for servicing purposes, use the Service menu to directly active or deactivate included and optional valves and vacuum pump.

Valves

In this menu the connected valves can be opened or closed individually.


Pump

User can start or stop the pump manually without actual distilling.

Temperature adjustment

User can set the different temperature to zero when there is no vapor liquefaction. The correct temperature adjustment helps the auto distillation mode to start and stop the distilling.

8. Settings

Languages

The "Languages" option allows the user to select the desired language.

Units

The "Units" option allows the user to select the desired unit for displaying the pressure and the amount.

Display

The "Display" option allows the user to change the background color and brightness of the working screen, and display "Firmware Update Info".

Sound

The "Sound" option allows the user to set the volume and activate/deactivate the key tone.

Factory Settings

This option allows the user to reset Menu values, Programs, Solvent library, or all the system settings to the original values at time of delivery.

Communication

Device name: The device name is used to identify a device.

Bluetooth: The "Bluetooth" option allows the user to activate/deactivate the "Bluetooth" function.

Labworldsoft 5 protocol: To use the RV 10 auto in Labworldsoft 5 this option has to be active.

Information

This option offers the user an overview of the most important system settings of the device.

Setting the lower end stop

CAUTION

Depending on the size of the piston, the setting angle of the rotation drive and the position of the heating bath and lift, the evaporation flask can be in contact with the bottom of the heating bath.

CAUTION

Glass may break!

Limit the lower lift position with the variable end stop.


The lift can be moved to any desired position in manual mode using the “▲” and “▼” keys. There is no automatic cut-off if a collision occurs.

- Press the “▼” key until the lift reaches the desired position.

Note: The evaporation flask should be 2/3 immersed in the heating bath.

- To adjust the stop element (1), press the center button (2) on the front of the lift.
- Move the stop (1) to the desired position
- Press the “▲” key until the lift has reached the upper end stop.

Note: The path is limited from 0 – 6 cm.

Check that the limit stop is functioning properly:

- Lower the lift by holding down the “▼” key.
- The drive stops automatically when the preset desired lower end position has been reached.
- Move the drive back to the upper position.

To activate the safety lift following a long period of inactivity, use the motor to move the lift to the top or bottom position several times prior to starting the distillation. (See the section entitled Safety Instructions - Safety Lift!)

Fill the evaporation flask

Manual filling: You can fill the evaporation flask manually prior to creating the vacuum. The evaporation flask should not be filled more than half its volume.

Automatic filling: Prior to filling the evaporation flask, a vacuum controller is used to regulate the glass apparatus to the target pressure.

- Now fill the evaporation flask using the back feed line.
- Due to the vacuum present, the solvent is suctioned into the evaporation flask. This enables you to keep solvent loss due to suctioning to a minimum.

CAUTION

The maximum allowed load (evaporation flask and contents) is 3 kg.

Setting up the heating bath

CAUTION

Refer to the operating instructions for the heating bath IKA HB!

- Move the lift to the bottom position and check the position of the heating bath in relation to the evaporation flask. When using larger evaporation flasks (2 or 3 litres) or depending on the angle of the rotation drive, you can move the heating bath 50 mm to the right.
- Fill the heating bath with the tempering medium until the evaporation flask is surrounded by tempering medium to 2/3 of its volume.
- Switch on the rotation drive and slowly increase the speed

Note: Avoid creating waves.

- Switch on the heating bath using the main on/off switch.

Note: Avoid stress on the glass due to different evaporation flask and heating bath temperatures when lowering the evaporation flask into the heating bath!

Note: If non-original accessories are used that are not supplied by IKA, then it is possible that the 50 mm travel range provided by the heating bath will not be sufficient. This applies in particular when using 3 L evaporator flasks with foam brake. Use the IKA RV 10.3000 extension plate to extend the heating bath travel range by 150 mm.

Interfaces and outputs

The device can be operated in "Remote" mode via an RS 232 or USB interface using *labworldsoft®* laboratory software.

The RS 232 interface at the back of the device, is fitted with a 9-pole SUB-D jack and can be connected to a PC. The pins have serial signals. The USB interface is located on the left-hand side of the display on the drive unit and can be connected to a PC using the USB cable supplied.

Note: Please comply with the system requirements together with the operating instructions and help section included with the software.

USB Interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized. Use the USB interface in conjunction with *labworldsoft®* for operation in "Remote" mode and also to update the firmware.

Installation:

First, download the latest driver for **IKA** devices with USB interface from: <http://www.ika.com/ika/lws/download/usb-driver.zip>. Install the driver by running the setup file. Then connect the **IKA** device through the USB data cable to the PC. The data communication is via a virtual COM port. Configuration, command syntax and commands of the virtual COM ports are as described in RS 232 interface.

Serial interface RS 232

Configuration:

- The functions of the interface connections between the stirrer machine and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s.
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request.

Command syntax and format:

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the device (Slave).
- The lab device sends only at the computer's request. Even fault indications cannot be sent spontaneously from the machine to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with CR LF (Code: hex 0x0d hex 0x0A) and have a maximum length of 128 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific **IKA** commands serve only as low level commands for communication between the rotary evaporator and the PC. With a suitable terminal or communication program these commands can be transmitted directly to the rotary evaporator. The **IKA** software package, *labworldsoft®*, provides a convenient tool for controlling the rotary evaporator and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarizes the (NAMUR) commands understood by the **IKA** control equipment.


Abbreviations used:

- m = Numbering parameter (integer)
- X = 2 Temperature heating bath
- X = 3 Safety temperature heating bath
- X = 4 Rotation speed
- X = 60 Interval time (1 - 5999 seconds, 1 <= m >=5999)
- X = 61 Timer (1 - 1440 minutes, 1 <= m >=1440)
- X = 62 Upper lift position (OUT_SP_62 1-> drive lift up)
- X = 63 Lower lift position (OUT_SP_62 1-> drive lift down)
- X = 66 Vacuum controller pressure
- X = 70 Vacuum controller hysteresis
- X = 74 Tempering medium heating bath (OUT_SP_74
0=oil, OUT_SP_74 1=water)

NAMUR Commands	Function
IN_NAME	Request designation
IN_PV_X X = 4	Read actual value
IN_SOFTWARE	Request software Id-number, date and version
IN_SP_X X = 4	Read target value input
OUT_SP_X m X = 1,60,61,62	Set target value to m
RESET	Switch to normal operation
START_X X = 4,60,61,62	Switch on appliance (remote) function
STATUS	Status output 0:Manual operation without interruption 1:Automatic operation Start (without interruption) ERROR z (z error number see table)
STOP_X X = 1,60, 61, 62	Switch off appliance function. The variable set with OUT_SP_X remains


PC 1.1 cable

This cable is required to connect the 9-pin connector to a PC.

USB 2.0 cable

This cable is required to connect USB interface to a PC.

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

The seal on the glass condenser should be checked at regular intervals and replaced if necessary.

Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by IKA to clean IKA devices.
These are water (with tenside) and isopropanol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com,
- software version

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Note: For glassware unit, please contact your local dealer and do not send the glassware back to our factory.

Safety lift

The safety lift must be checked daily prior to operation!

After a long period of inactivity (approx. four weeks), use the motor to move the lift to the top or bottom position several times prior to starting distillation. (See the section entitled Safety Instructions - Safety Lift)! If the safety lift is not working, please contact the **IKA** Service department.

Accessories

For further accessories see www.ika.com.

Error messages

Any malfunctions during operation will be identified by an error message on the display.

Once a serious error message has been displayed, the lift moves to the top end position and the device can no longer be operated.

Proceed as follows in such cases:

- switch off device using the mains switch,
- carry out corrective measures,
- restart device.

Error message	Effect	Cause	Correction
No pressure change	Pressure gradient analysis does not show any deviation	Pump is not running Vacuum hose isn't connected RV 10.4002 valve control cable is disconnected Failure in hose system Recipient isn't closed	Check pump control cable connection Check pump power supply Switch on pump Check vacuum hose connection (pump, recipient) Check RV 10.4002 valve control cable connection Check correct vacuum hose connection (IN / OUT/ VENT) Close recipient
System not tight	Pressure gradient analysis shows pressure deviation, but set value cannot be reached. Pressure set value cannot be reached.	Vacuum hose connection isn't tight Recipient isn't tight Vacuum pump power isn't sufficient Venting valve leakage	Check vacuum hose connection Check recipient Check technical data of the pump Contact service department
Pressure out of range	Vacuum is out of range Vacuum is bigger than room pressure	Pressure in the recipient is too high Sensor is broken.	Check the air flow of pump and vent recipient Contact service department
Vacuum sensor error	Vacuum sensor output is too low	Sensor is not connected Sensor is broken.	Check the sensor connection Contact service department
Vacuum calibration error	Vacuum calibration value is out of range	Input wrong vacuum actual value when doing calibration. Set wrong vacuum calibration value Sensor is broken	Check and calibration again Check the set calibration value Contact service department
Boiling point detecting error	Boiling point cannot be detected. Differ-temperature analysis does not show boiling point	Heating bath does not heat up Solvent with extreme low boiling point is used	Check heating bath Solvent distillation only manually
Venting error	Pressure gradient analysis does not show any deviation after pressing the "Venting" button.	Recipient wasn't evacuated Failure in hose system Venting valve doesn't work Front foil button is damaged	Evacuated recipient Check correct vacuum hose connection (IN / OUT/ VENT) Contact service department Contact service department
Analog pump not connected	Analog pump is disconnected		Check the analog pump cable connection
No rotation	Rotary drive doesn't start or exhibits very large speed variation	Rotary drive is overload or blocked Motor or motor cable is damaged	Reduce the load Contact service department
Temperature sensor not connected	Differential temperature shows big minus value	Temperature sensor is not connected Temperature sensor is damaged	Check the Temperature sensor connection Contact service department
No temperature difference	No differential temperature increase of cooling water	No solvent or too less solvent is distilled No water flow Water flow is wrongly connected	Check the solvent Check flow of cooling water Check flow direction of cooling water
Temperature out of range	Differential temperature out of the measurement range	Too fast distilling Temperature sensor is damaged	Increase the set vacuum Contact service department

Error message	Effect	Cause	Correction
Temperature calibration error	Temperature calibration value is out of range	Wrong simulator resistance is chosen PCB is broken	Check the simulator resistance Contact service department
Temperature adjustment error	Water in and water out sensor temperature have big deviation	The cooling water temperature is changing Temperature sensor is not calibrated Temperature sensor is broken	Check the chiller and wait until the chiller reaches the set temperature Make sure there is no distilling Contact service department Contact service department
Inside temperature too high	Analysis of internal PCB temp. Sensor. Inside temperature reaches limit value.	Room temperature > 40°C PCB is broken	Switch off device to cool down Contact service department
Flow rate out of range	Water flow rate of cooling water is out of range	Water flow rate is too high Water flow sensor is damaged	Reduce the water flow rate Contact service department
Flow rate too low	Water flow rate of cooling water is less than the minimum limit value	Water flow rate is too low Water flow sensor is damaged	Increase the water flow rate Contact service department
Condenser overload	Cooling power limit value is exceeded	Condenser limit value is too low Distilling is too fast	Increase the condenser overload limit Increase the set vacuum
Distillation dry	There is no more solvent distilled	Distillation is completed	Stop the distilling
Lift error	Lift can not reach the end position	Lift is overloaded or lift is blocked Position sensor failure Lift motor, cable or PCB is damaged	Check the lift Contact service department Contact service department
Internal communication error	Internal Communication Error	Internal communication error	Switch device on/off, Contact service department
PC communication error	Error command received from PC via RS 232 or USB interface	LWS uses wrong device for RV 10 control Wrong tools are used to communicate with RV 10 control Unstable connection	Check the LWS setting Check the PC tools Check the connection
Heating bath communication error	Communication with heating bath is not available	Infrared interface is blocked Heating bath is switched off or detect error Heating bath with infrared interface is not present PCB is damaged	Check and clean the infrared interface Switch on heating bath Check the heating bath Contact service department
EEPROM error	Error while reading or writing the EEPROM	EEPROM is defective	Contact service department
Storage error	Read or write internal flash error	Read or write internal flash error	Contact service department

If the actions described fail to resolve the fault or another error message is displayed, then take one of the following steps:

- contact **IKA** service department,
- send the device for repair, including a short description of the fault.

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine directly to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Note: For claims on glassware units, please contact your local dealer and do not send the glassware back to our factory.

Technical Data

Operating voltage range	Vac	(100...240) ± 10%
Rated voltage	Vac	100...240
Frequency	Hz	50/60
Power input without heating bath	W	100
Power input operation "standby"	W	3.3
Speed	rpm	0/5...300
Speed tolerance	rpm	± 1 (Set speed < 100 rpm)
	%	± 1 (Set speed ≥ 100 rpm)
Speed display		digital
Dimensions of visible display area (W x H)	mm	70 x 52
Display		TFT-Display
Multiple languages		yes
Right and left movement/ interval operation		yes
Smooth start		yes
Lift		automatic
Lifting speed	mm/s	50
Stroke	mm	140
Setting of lower end stop	mm	60, contactless
Head angle adjustable		0°...45°
Timer	hh:mm:ss	00:00:00...99:59:59
Interval	mm:ss	00:00...60:00
Integrated vacuum controller		yes
Measurement range vacuum	mbar (hPa)	1...1050
Adjustment range vacuum	mbar (hPa)	2...1014
Measurement accuracy vacuum	mbar	± 2 (± 1 Digit) if calibration is done at a constant temperature
Adjustment accuracy vacuum	mbar	adjustable hysteresis (Two-position vacuum control)
	mbar	1 (Speed-vacuum control, set value < 100 mbar)
	%	1 (Speed-vacuum control, set value ≥ 100 mbar)
Measurement range of differential temperature	K	7
Offset balancing range	K	± 0.5
Interface		USB, RS 232
Ramp programming		yes
Distillation processes programmable		yes
Remote control		with accessory <i>labworldsoft</i> ®
Cooling surface (Standard condenser RV 10.1, 10.10, 10.2, 10.20)	cm²	1500
Minimum cooling water flow rate	l/h	30
Maximum cooling water flow rate	l/h	100
Cooling water pressure	bar	1
Cooling water temperature range	°C	18...22, constantly
Fuse		2 x T1, 6 A 250 V 5 x 20
Perm. On-time	%	100
Perm. ambient temperature	°C	5...40
Perm. humidity	%	80
Protection acc. to DIN EN 60529		IP 20
Protection class		I
Overvoltage category		II
Contamination level		2
Weight (no glassware; no heating bath)	kg	18.6
Dimensions (W x D x H)	mm	500 x 440 x 430
Operation at a terrestrial altitude	m	max. 2000 above sea level

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited