IKA

designed for scientists

IKA LR-2.ST

BETRIEBSANLEITUNG	DE	3
OPERATING INSTRUCTIONS	EN	14
MODE D'EMPLOI	FR	25
INDICACIONES DE SEGURIDAD	ES	38
VEILIGHEIDSinSTRUCTIES	NL	40
NORMI DI SICUREZZA	IT	42
SÄKERHETSANVISNINGAR	SV	44
SIKKERHEDSHENVISNINGER	DA	46
SIKKERHEDSHENVISNINGER	NO	48
TURVALLISUUSOHJEET	FI	50
Υπόδειξεις ασφάλειας	EL	52
INSTRUÇÕES DE SEGURANÇA	PT	54

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	Page
Declaration of conformity	14
Explanation of symbols	14
Safety instructions	14
Correct use	15
Unpacking	16
Assembling	16
LR-2.ST stand system	17
Stirrer	18
Reactor head	18
Reaktor vessel and clamping mechanism	19
Safety disconnection LR-2.SI	19
Use of dispergiersing element	21
Maintenance and cleaning	21
Cleaning	21
Sterilization	22
Maintenance	22
Accessories	23
Materials in contact with medium / laboratory reactor	24
Guarantee	24
Technical data	24

Declaration of Conformity

We declare under our sole responsibility that this product corresponds to the regulations 2006/42/EC and conforms with the standards or standardized documents EN ISO 12100.

Explanation of symbols

 DANGER	Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
 WARNING	Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
 CAUTION	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 NOTICE	Indicates practices which, if not avoided, can result in equipment damage.
 CAUTION	Danger! - note on hazards arising from a hot surface.
 CAUTION	Indicates crushing risk of fingers/hand.

Safety instructions

General instructions

Read the operating instructions completely before starting up and follow the safety instructions.

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Implement the accident prevention regulations for this type of activity and workplace. Wear your personal protective equipment at all times when working with the laboratory reactor.
- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Do not use flammable or explosive gases, vapours or solvents in the LR-2.ST laboratory reactor!
- Any reaction vapour released at the reactor lid must be extracted.
- The appliance is not designed for use in explosion hazard (Ex-certificate) areas!
- Certain applications and materials may be hazardous. You should take precautions to prevent contact with, or inhalation of, toxic liquids, gases, fumes, vapours or powders.
- Risks may also be posed by biological or microbiological substances.

DANGER

Do not touch rotating parts. Moving and rotating equipment parts also constitute a hazard.

CAUTION

There is a danger of crushing when lowering the reactor lid.

Power supply/switching off the device

CAUTION

The safety disconnecter does not apply to the heating bath.

- Appliances that need to continue operation once the reactor cover has been opened (e.g. thermostat), must not be operated via the safety disconnecter.
- Only operate the device under a laboratory extractor hood (closed on all sides) or other appropriate safety device.
- There must be an EMERGENCY STOP facility which disconnects all poles to the laboratory reactor system from the power supply. The laboratory reactor must only be operated using the safety disconnecter provided or in an extractor hood with emergency stop.
- Connect the laboratory stirrer and the disperser drive to the mains power supply over the LR-2.SI safety disconnecter. Connect these appliances to the multisocket on the safety disconnecter only
- Connect these appliances to the locked socket strip with safety disconnecter only.

Performing trials

DANGER

Chilblains/burns due to low/high temperatures on the hose coupling. Open the hose coupling only at room temperature.

- The laboratory stirrer/disperser is fitted with a limit switch which prevents the unit from being started until it is in the fully lowered position. The anchor stirrer and the dispersing element must only be operated when the reactor vessel is fully closed.
- Do not start up the laboratory reactor when it is open.

Device configuration

DANGER

Glass breakage due to pressure differences at the quick-release connectors. When using quick-release connectors the double jacket is completely sealed. Pressure differences can occur due to cooling/warming of the medium.

WARNING

Tempering liquids must be at room temperature when changing/releasing. The equipment is not designed for overpressure use.

- The reactor system must always be ventilated when working under normal pressure in order to prevent any pressure build-up caused by highly volatile gases or unpredictable reaction pressure gradients. Condense volatile gases using a cooler with a standard ground connector (e.g. a reflux condenser) on the reactor cover.
- The appliance is designed for operation at vacuum settings up to 25 mbar (see Accessories).
- Please follow the safe procedures specified for handling reactants and reaction products.
- The LR 2000.1 and LR 2000.2 laboratory reactor vessels may only be heated without pressure using the double shell with thermostats or similar heat sources. Do not under any circumstances use a heating mantle or a hot plate (Bursting hazard!).
- We recommend that you use the **IKA®** HBR 4 digital heating bath to temper the LR 2.1 single-walled laboratory reactor vessel.

Heating

CAUTION

Burns Hazard. The heated reactor vessel and cover are very hot. Wear protective gloves when handling hot parts.

- Ensure that the thermostat used for tempering is fully efficient. A defective thermostat can cause uncontrolled reactions.

- The maximum permitted temperature inside the reactor vessel is 230 °C. Do not exceed this temperature.
- Before you fill the reactor vessel, ensure that the reagents used do not corrode the seal!

Stand

Undo the screws on the stand carefully. Without counterbalance the stand can move upwards uncontrollably.

Accessories

- Only use **IKA**® approved accessories!
- Use only original **IKA**® spare parts!
- Protect the instrument and accessories from bumping and impacting.
- Please observe the permitted speed for the stirring element. Never set higher speed.

Assembly

Ensure that the stirring tool is securely fixed.

Broken glass! Fit the hose coupling carefully.

- Fix all hoses in the laboratory extractor hood.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
 - stress due to incorrect assembly,
 - external mechanical hazards,
 - local temperature peaks.

Protection of the device

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact) supply.
- Covering or parts that are capable of being removed from the unit without accessory equipment have to be reattached to the unit for safe operation in order to prevent, for example, the ingress of fluids, foreign matter, etc..
- The appliance may only be opened by experts.

Correct use

NOTICE

The anchor stirrer must only be operated at a speed of max. 300 rpm.

You can use these in conjunction with labworldsoft® laboratory software to create a high-performance system for measuring, controlling, regulating and fully documenting all laboratory processes.

Sample applications:

- Manufacture of cremes, lotions, emulsions and liposome preparations in the pharmaceutical and cosmetic sector (lipstick, sun cream, eyeliner).

- Mixing of solids such as calcium carbonate, talc and titanium oxide with liquid polymers.
- Mixing of additives and solid polymer compounds with mineral oils.
- Dispersion and mixing of solids and fibres in liquids and polymers.
- Manufacture of milk products and dietary products.
- Extraction of enzymes from biomass.

Unpacking

Unpacking

- Unpack the device carefully
- Any damage should be notified immediately to the shipping agent (post office, railway network or transport company)

The LR-2.ST laboratory reactor system is supplied as a fully functional unit in various configurations depending on customer requirements. Accessories, such as ULTRA-TURRAX, temperature gauges and reactor vessels are supplied separately.

Scope of delivery

- LR-2.ST stand system complete with reactor cover
- Overhead stirrer
- LR-2.SI safety disconnecter

Assembling

LR-2.ST stand system

NOTICE

The stand system holds the drive unit complete with the reactor cover and the stirring instrument, the reactor vessel and the safety disconnecter.

An adjustment screw on the rear, right-hand side foot of the stand can be used to level the stand and compensate for any unevenness in the mounting surface.

Undo the lock nut. Turn the knurled screw clockwise to slightly raise the foot of the appliance, turn it anticlockwise to lower it. When the appliance is level, tighten the lock nut.

M8 cheese-head screw with baseplate and side support

Adjuster screw for levelling

If you are working with the LR 2.1 single-walled reactor vessel and the HBR 4 digital heating bath, you must remove the baseplate including clamping bracket and one side support from the stand.

Undo the top two M8 cheese-head screws using the DIN 911-6 cranked wrench key provided. Remove the baseplate, undo the bottom cheese-head screws and remove the side support.



The remaining side support is used to fit the safety disconnecter. Place the plastic cover onto the aluminium profile.

Undo the knurled screw to adjust the height

To raise the whole assembly by hand, first undo the two knurled screws on the left and right of the stand legs. The pressure springs built into the telescopic stand mean that the whole reactor can be easily raised and lowered. The laboratory stirrer mount with integrated stirrer shaft bearing is fitted to the stand adapter (connected to the stand legs) with two tie rods. The ball-bearing, vacuum-tight stirrer shaft is connected to the stirrer shaft drive by a coupling flange. The stirrer shaft features covered

precision bearings and is sealed with rotary shaft seals with PTFE sealing lips for vacuum operation. Penetration of the stirrer shaft inevitably produces wear debris from the stirrer shaft seal. In order to keep this debris away from the product, a protective ring has been fitted to the stirrer shaft from below. The anchor stirrer is also inserted into the stirrer shaft from below. The torque is transferred by a square shaft to the anchor stirrer. This is secured axially by means of a retaining ring.

Stirrer

NOTICE

Refer to the operating instructions for the laboratory stirrer. See "Technical data" section

The laboratory stirrer is a fixed component of the laboratory reactor system.

The following IKA® laboratory stirrers are available for this system:

- EUROSTAR 200 control P4
- EUROSTAR 200 control
- EUROSTAR 100 control
- EUROSTAR 60 control

Reactor head

DANGER

**Hazardous pressure build-up!
Ensure that the reactor vessel is properly vented and can release sufficient air!**

The reactor cover is screwed onto the stirrer shaft bearing using a thread. The reactor cover can be easily removed for cleaning or sterilisation purposes (see Cleaning / Sterilisation).

An O-ring provides the seal between the stirrer shaft bearing and the reactor cover.

The reactor cover is fitted with three NS 29/32 standard ground connectors and two NS 14/23 standard ground connectors. These connectors are screwed to the cover and sealed by O-rings.

One of the NS 14/23 connectors is fitted with a ventilation and air vent valve, whilst the other connectors are closed off with DIN 12252 ball valves.

All sealing rings on the reactor cover are made of FPM (perfluoroelastomer).

Dispersing elements (ULTRA TURRAX), temperature sensors, flow breakers and other accessories can be fitted to the free connectors. In order to do this, unscrew the threaded connector plug from the cover using the pin provided and then fit the device (e.g. LR 2000.60 temperature sensor) in its place.

The lock ring used to secure the reactor is located on the rim of the cover.

The clamping force of the lock ring can be set when it is open using the adjuster screw.

Adjuster screw

Reactor vessel and clamping mechanism

Configuration with single-walled reactor vessel

NOTICE

Secure the reactor vessel using the clamping strap on the stand. Screw the hose adapters to the tempering hoses.

The single-walled reactor vessel is directly centred with and connected to the reactor cover using the lock ring.

Position the stainless steel cover plate on the tempering bath and centre the vessel with the cork ring underneath the reactor cover. Lower the stand as far as the top glass edge and secure the glass vessel by closing the clamping lever on the lock ring.


When you insert the vessel, check that the sealing contact surface of the glass vessel is undamaged and that the seal on the reactor cover fits properly.

Configuration with double-walled reactor vessel

The LR 2000.1 double-walled reactor vessels without bottom discharge valve and LR 2000.2 double-walled reactor vessels with bottom discharge valve are set onto the cork ring provided.

1. Connect the metal adapter on the glass tubing connector to the coupling (quick-release) LT 5.25.

Quick-release connector

Metal adapter on glass tubing connector

2. Connect the LT 5.25 to the tempering hoses.
3. Then connected the fully assembled hoses to the circulation thermostat.

Lower the drive unit including reactor cover and anchor stirrer and connect the reactor vessel to the cover by closing the lock ring. Firstly check that the sealing surface on the glass vessel is intact and that the seal on the reactor cover fits properly.



Clamping strap

Double-walled vessel

NOTICE

The IKA® LR-2.SI safety disconnecter disconnects the laboratory reactor system from the mains power supply when the red EMERGENCY STOP button is pressed. The safety disconnecter must be actively switched on again if the power supply is interrupted. The appliances cannot automatically restart.

You can mount the safety disconnector power box on either the left or the right aluminium profile using two slider blocks. Insert the slider blocks into the profile groove with the spring-loaded ball to the back. You do not need to remove the baseplate. Position the slider blocks in the profile in such a way that you can screw on the power box with the two M6 cheese-head screws using the SW5 Allen key.

Slide the wireless safety switch down to the bottom end of the top (thicker) stand rod and secure it using the clamping screw.

The assembled safety switch should be pointing downwards. Turn the bracket so that it is positioned in such a way that it will not prevent you from operating the clamping screw for the stand properly.

Lower the stand until the reactor cover is sitting on top of the reactor vessel and secure it in this position by tight-

ening the knurled screws.

Fit the second safety switch (with cable) to the bottom (thinner) stand rod in this position and secure it so that it is approx. 5 mm away from the top safety switch. Align the two switches with each other.

Attach the socket bracket to the top of the stand and place the multi-socket on top.

Connect the mains plugs of the stirrer drive and the disperser drive (if you are using one) to the multi-socket on the safety disconnector.

Fix the power cables together with cable ties.

Connect the power supply cable of the safety disconnector to your power supply.

Switch on the laboratory stirrer drive at its main switch. Close the reactor cover and press the green button on the safety box. The stirrer drive will start up.

To switch off the stirrer drive, press the red button or open the reactor cover.

The limit switch trigger on the stand rod (switch distance approx. 5 mm) also breaks the circuit if the reactor cover is open. This prevents any danger from rotating elements such as anchor stirrers or dispersing elements.

The safety disconnector must be switched back on again each time it is switched off (use the green button). After a power failure it will also be necessary to switch on the safety disconnector.

Use of dispersing element

NOTICE

Use of a dispersing element is only permitted when the reactor contains a working volume of at least 800 ml. Using a dispersing element with volumes lower than this will result in dry running which will damage the axial face seal in the dispersing element.

The following dispersing elements can be used with the ULTRA TURRAX T 25 drive unit:


All parts in contact with the medium are made of 1.4571 stainless steel or FPM/SIC. The minimum/maximum immersion depth is 40/225 mm. These dispersing elements can be used in the reactor system at vacuums of up to 25 mbar.

The LR 2000.40 shaft holder is necessary to install the dispersing element in the reactor cover.

	S 25 KV - 18 G	S 25 KV - 25 G	S 25 KV - 25F	S 25 KV 2802
Rotor/stator gap width (mm)	0,3	0,5	0,5	-
Max. temperature (°C)	220	220	220	220
Suspension fineness (µm)	10-50	15-50	5-25	-
Emulsion fineness (µm)	1-10	1-10	1-5	-

Assembly:

- Insert the circlip into groove on the shaft
- Place shaft into the shaft holder from below
- Insert the circlip into upper groove
- Tighten clamping screw

The circlips required are supplied with the laboratory reactor.

The assembled shaft holder is replaced by a standard ground connector in the reactor cover.

Commissioning

NOTICE Follow the operating instructions given in this manual!

You should also refer to the operating instruction manual provided with the other items of equipment delivered to you!

The reactor system must be operated under a closed fume hood or behind a protective screen!


The reactor must only be used with the safety disconnecter supplied. The safety disconnecter ensures that the system can be completely disconnected from the mains power supply when required!

Only connect the laboratory stirrer and the disperser drive to the multi-socket. The multi-socket is mounted on the back of the equipment and is connected to the safety disconnecter.

Fit the safety disconnecter!

Any reaction vapour released at the reactor cover must be extracted!

The appliance will not restart automatically after a power failure. You must switch the safety disconnecter back on.

 The socket used must be earthed (fitted with earth contact).

- Secure the reactor vessel using the lock ring on the reactor cover.
- Connect the power supply cables to the mains supply.
- Switch on the safety disconnecter by pressing the green button.
- Start the laboratory stirrer.


The **IKA® LR-2.SI** has an EMERGENCY STOP button and also a safety disconnecter which provides an alternative means of shutting down the power supply circuit when the reactor cover is open. This protects you from rotating elements such as anchor stirrers or dispersing elements.

The glass reactor vessel must always be fully immersed in the tempering unit. With this equipment configuration you must ensure that any vapours released from the tempering bath are extracted.

Maintenance and cleaning

Cleaning

NOTICE When the laboratory reactor is open the reactor vessel can be completely removed from the reactor so that the interior can be cleaned.

 Remove the device from the mains before cleaning.


The axially inserted anchor stirrer can be removed from the stir shaft with a quick downward pull. The Teflon scrapers (LR 2000.10 anchor stirrer) can be removed from the anchor stirrer for thorough cleaning.

If no additional devices such as ULTRA TURRAX, temperature sensors or flow breakers are installed, the reactor cover can be unscrewed from the stirrer shaft bearing for cleaning.

First remove the connector components, the connectors for the additional devices and pull the protective ring off the stir shaft on the cover. The stirrer shaft and the protective ring can also be cleaned now.

All O-rings can now be carefully removed by hand or using a blunt tool.

Removing the anchor stirrer

Sterilization

Laboratory reactor components may only be sterilised using these methods

Chemical methods:

Germicidal solutions (formalin, phenol, alcohol etc.) are good general-purpose disinfectants. It is important that all disinfectant residues are then removed with sterilized water.

Sterilization with moist heat:

This denotes the use of a steam jet pressurized to 2 bar at 120 °C.

Sterilization with hot air:

Hot-air sterilization is usually performed at 160 °C or 190 °C.

Maintenance

Refer to the operating instructions for individual components.

The IKA® laboratory reactor system requires little maintenance. The ball bearings are sealed and lubricated for life.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name and item-no. of the spare part.

Before you send any equipment to us for repair, ensure that it has been cleaned of any toxic or corrosive substances. Use only water with a detergent additive containing a surface active agent or, in cases of extreme contamination, isopropyl alcohol for cleaning parts not in contact with the medium.

Repairs

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For this, use the “**certificate of compliance**” form which you can obtain from IKA® or can download a version for printing from the IKA® website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Materials in contact with medium/laboratory reactor

Reactor cover	Stainless steel 1.4571	Scraper	PEEK
Reactor vessel	Borosilicate glass 3.3	O-ring	FFPM
Anchor stirrer	Stainless steel 1.4571	Shaft seal	PTFE/1.4571

Guarantee

In accordance with **IKA®** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual

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We declare under our sole responsibility that this product corresponds to the regulations 2006/42/CE and conforms with the standards or standardized documents EN ISO 12100.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.