

IKA

designed for scientists

IKA MIDI MR 1 digital IKA MAXI MR 1 digital

| | | | | | |
|---------------------------------------|----|----|-----------------------------|----|-----|
| Betriebsanleitung Ursprungssprache | DE | 4 | | | |
| Operating instructions | EN | 11 | Használati utasítás | HU | 103 |
| Mode d'emploi | FR | 18 | Navodilo za delovanje | SL | 110 |
| Руководство по эксплуатации | RU | 25 | Návod na použitie | SK | 117 |
| Instrucciones de manejo | ES | 33 | Kasutusjuhend | ET | 124 |
| Handleiding | NL | 40 | Lietošanas instrukcija | LV | 131 |
| Istruzioni per l'uso | IT | 47 | Darbo instrukcija | LT | 138 |
| Driftsanvisning | SV | 54 | Ръководство за експлоатация | BG | 145 |
| Driftsinstruks | DA | 61 | Instrucțiuni de utilizare | RO | 152 |
| Driftsvejledning | NO | 68 | Οδηγίες χρήσης | EL | 159 |
| Käyttöohje | FI | 75 | 使用说明 | ZH | 166 |
| Instruções de serviço | PT | 82 | 取扱説明書 | JA | 172 |
| Instrucja obsługi | PL | 89 | 사용 설명서 | KO | 179 |
| Návod k provozu | CS | 96 | | | |

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Contents

| | Page |
|--------------------------------|------|
| EU Declaration of conformity | 11 |
| Explication of warning symbols | 11 |
| Safety instructions | 11 |
| Correct use | 13 |
| Unpacking | 13 |
| Commissioning | 13 |
| Interface and output | 14 |
| Maintenance and cleaning | 16 |
| Error codes | 16 |
| Accessories | 16 |
| Technical data | 17 |
| Warranty | 17 |

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at sales@ika.com.

Explication of warning symbols

| | |
|--|--|
|  DANGER | Indicates an (extreme) hazardous situation, which, if not avoided, will result in death, serious injury. |
|  WARNING | Indicates a hazardous situation, which, if not avoided, can result in death, serious injury. |
|  CAUTION | Indicates a potentially hazardous situation, which, if not avoided, can result in injury. |
|  NOTICE | Indicates practices which, if not avoided, can result in equipment damage. |
|  NOTICE | Note the hazards of magnetism! |

Safety instructions

General information:

- **Read the operating instructions completely before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Socket must be earthed (protective ground contact).
- The socket for the mains cord must be easily accessible.
- The feet of the device must be clean and undamaged.
- Check the device and accessories for damage before each use them. Do not use damaged components.
- Do not use the device in explosive atmospheres, it is not EX-protected.
- With substances capable of forming an explosive mixture,

appropriate safety measures must be applied, e.g. working under a fume hood.

- To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.
- The device may heat up when in use.
- Do not move or transport the device when it is operating or connected to the power supply.
- Safe operation is guaranteed only with the use of original IKA accessories.
- Always disconnect the plug before attaching accessories.
- The device can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The device can only be opened by technical experts, even during repair. The device must be unplugged from the power supply before opening. Live parts inside the device may still be live for some time after unplugging from the power supply.

NOTICE

Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).

- The working plate can heat up due to the action of the drive magnets at high motor speeds.
- The centre of gravity of the assembly must lie within the surface on which it is set up.
- It may be possible for wear debris from rotating accessory parts to reach the material being processed.
- When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

(Source: *Römpps Chemie-Lexikon* and *"Ulmann", Volume 19*)

For protection of the user:

WARNING

Wear your personal protective equipment in accordance with the hazard category of the media

- to be processed. There may be a risk from:
- splashing and evaporation of liquids
 - ejection of parts
 - release of toxic or combustible gases.

DANGER

Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.

DANGER

Process pathogenic materials only in closed vessels under a suitable fume hood. Please contact **IKA** if

you have any questions.

WARNING

Beware of hazards due to:

- flammable materials
- incorrect container size
- overfilling of media
- unsafe condition of container.

For protection of the device and accessories:

- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The voltage stated on the type plate must correspond to the mains voltage.
- Removable parts must be refitted to the device to prevent the infiltration of foreign objects, liquids etc..
- Protect the device and accessories from bump and impact.
- Do not cover the device, even partially e.g. with metallic plates or film. This may result in overheating.

Performing trials:

NOTICE

Reduce the speed if:

- the medium splashes out of the vessel
- the device is not running smoothly
- dynamic forces start to cause the device and/or the vessels placed on it to move around.

- Make certain that the device is set at the lowest speed (left hand position) before commissioning; otherwise, the device will start to run at the speed set in last operation. Gradually increase the speed.
- After an interruption in the power supply or a mechanical interruption during working process, the device will restart automatically.

Correct use

The **IKA MIDI MR 1 digital** and **MAXI MR 1 digital** is a magnetic stirrer without heating function. The devices are suitable for stirring and mixing liquid substances up to 50 liter (**MIDI MR1 digital**) and 150 liters (**MAXI MR1 digital**). The built-in stirring drive permits the stirring of substances with the aid of a magnetic stirring bar in the vessel. The mixing intensity is dependent on the motor speed and the size of the magnetic bar.

• Use:

- for stirring and mixing liquids.

Intend use: Tabletop device

• Area of use:

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

The safety of the user cannot be guaranteed:

- if the device is operated with accessories that are not supplied or recommended by the manufacturer,
- if the device is operated improperly or contrary to the manufacture's specifications,
- if the device or the printed circuit board are modified by third parties.

Unpacking

• Unpacking:


- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder).

• Delivery scope:

- **IKA MIDI MR 1 digital** or **MAXI MR 1 digital** according to order
- Magnetic stirring bar
- Mains cable
- USB 2.0 cable A – B
- Operating instructions
- Warranty Card.

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available.

The power socket used must be earthed (protective earthing).

If above conditions are met, the device is ready for operation after plugging in the mains plug.

If these conditions are not met, safe operation is not guaranteed and the device could be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under **"Technical Data"**.

The device can be operated in two modes, **"Continuous mode"** or **"Timer mode"**:


Fig. 3

Continuous mode:

- ☞ To switch the device on, turn the **"On/Off"** and **"Timer"** rotary knob (**B**, see **Fig. 3**) clockwise to the **"1"** position.
- ⇒ The device is now running in **"Continuous mode"**.
- ☞ To switch the device off, turn the **"On/Off"** and **"Timer"** rotary knob (**B**) to **"0"** position.

Timer mode:

- ☞ To switch the device on, turn the **"On/Off"** and **"Timer"** rotary knob (**B**) clockwise after the **"∞"** position.
- ⇒ The running time can be adjusted to any value between 0 and 55 minutes.
- ⇒ After the selected time period has elapsed, the switch will automatically return to the **"0"** position and the device will remain switched off.
- ⇒ The time period selected can be adjusted at any time.

Setting the motor speed:

- ☞ Set the motor speed by turning the **"Speed"** rotary knob (**C**, see **Fig. 3**). The motor speed can be set to between 0 to 1000 rpm for **MIDI MR 1 digital** and 0 to 600 for **MAXI MR1 digital**.

Interface and output

The device can be connected to a PC and operated with the laboratory software labworldsoft® through the RS 232 interface (G, see Fig. 2) or USB interface (H, see Fig. 2).

Note: Please observe the system requirements as well as the operating instruction and help section of the software.

USB interface:

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized.

Use the USB interface in conjunction with labworldsoft® for operation in "Remote" mode and also to update the firmware.

USB device drivers:

First, download the latest driver for IKA devices with USB interface from:

<http://www.ika.com/ika/lws/download/usb-driver.zip>.

Install the driver by running the setup file. Then connect the IKA device through the USB data cable to the PC.

The data communication is via a virtual COM port. Configuration, command syntax and commands of the virtual COM ports are as described in RS 232 interface.

RS 232 interface:

Configuration:

- The functions of the interface connections between the device and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.

- Transmission speed: 9600 bit/s.
- Data flow control: none.
- Access procedure: data transfer from the device to the computer takes place only at the computer's request.

Command syntax and format:

The following applies to the command set:

- Commands are generally sent from the computer (Master) to the device (Slave).
- The device sends only at the computer's request. Even fault indications cannot be sent spontaneously from the device to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific IKA commands serve only as low level commands for communication between the device and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the device. The IKA software package, Labworldsoft®, provides a convenient tool for controlling the device and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

Commands:

| Commands | Function |
|-------------|--|
| IN_NAME | Input description name. |
| IN_SP_4 | Reading the set rated value of rpm. |
| IN_PV_4 | Reading the real value of rpm. |
| OUT_SP_4 n | Setting the rated value of rpm to n. |
| OUT_SP_42@n | Setting the WD safety speed with the echo of the set value |
| OUT_WD1@m | Watchdog mode 1: When a WD1 event occurs, the stirring functions are shut down and message PC 1 is displayed. Set the watchdog time to m (20...1500) seconds, with echo of the watchdog time. This instruction starts the watchdog function and must be sent within the set watchdog time. |

| | |
|-----------|--|
| OUT_WD2@m | Watchdog mode 2: When a WD2 event occurs, the speed set point will be set to the WD safety set point speed. The PC 2 warning is displayed. The WD2 event can be reset with OUT_WD2@0-resetting also blocks the watchdog function. Set the watchdog time to m (20...1500) seconds, with echo of the watchdog time. This command starts the watchdog function and must be sent within the set watchdog time. |
| RESET | Switching off the device function. |
| START_4 | Starting the device (remote) function |
| STATUS_4 | Display of status: 10: Manual operation without fault 11: Automatic operation Start (without fault) 12: Automatic operation Stop (without fault) < 0: error code: - 1: error 1 - ... (see „Error codes“ table) |

Note: when disconnecting the RS 232 or USB cable while the device is running, the device stays in PC mode. To reset to function without PC, switch the device off and on again.

PC 1.1 Cable:

This cable is required to connect RS 232 interface (G) to a PC.


Fig. 4

USB 2.0 cable A - B:

This cable is required to connect USB interface (H) to a PC.


Fig. 5

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by IKA to clean IKA devices. These are water (with tenside) and isopropanol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.

- position number and description of spare part, see www.ika.com.

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

Any malfunctions during operation will be identified by an error message on the display.

Proceed as follows in such cases:

- ☞ Switch off device using the main switch at the back of the device.
- ☞ Carry out corrective measures.
- ☞ Restart device.

| Error code | Effect | Cause | Solution |
|---------------|--|-----------------------------|--|
| Err 10 | Stirring function doesn't start. | Speed setting signal error. | - Contact IKA service department. |
| Err 14 | Stop stirring or change according to relating setting. | PC communication failure | - Check communication cable. |

If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

Accessories

- **IKAFLON®** Magnetic stirring bar
- **TRIKA** Magnetic stirring bar
- **RSE** Stirring bar remover
- **PC 1.1** Cable
- **Labworldsoft®**

See more accessories on www.ika.com.

Technical data

| | | MIDI MR 1 digital | MAXI MR 1 digital |
|--|------------|--|-------------------|
| Operating voltage | VAC | 230 ± 10 % 115 ± 10 % 100 ± 10 % | |
| Frequency | Hz | 50 / 60 | |
| Power input | W | 70 | 80 |
| Power output | W | 19 | 35 |
| Max. load on top surface | kg | 75 | 200 |
| Stirred quantity max. (H ₂ O) | ltr | 50 | 150 |
| Maximum magnetic bar (L x Ø) | mm | 80 x 10 | 155 x 27 |
| Infinitely adjustable speed range | rpm | 0 ... 1000 | 0 ... 600 |
| Speed tolerance | | < ± 10 % of maximum speed | |
| Speed display | | LED | |
| Operating mode | | continuous operation and timer | |
| Timer | min | ∞ / 1 ... 55 | |
| Interface | | USB, RS 232 | |
| Fuse | A | 2 x T 4 A 250 V | |
| Permissible ambient temperature | °C | + 5 ... + 40 | |
| Permitted ambient humidity | % | 80 | |
| Permitted on-time | % | 100 | |
| Protection class according to EN 60 529 | | IP 21 | |
| Working plate dimensions (W x D) | mm | 350 x 350 | 500 x 500 |
| Working plate material | | stainless steel 1.4301 | |
| Dimensions (W x D x H) | mm | 360 x 430 x 110 | 505 x 585 x 110 |
| Weight | kg | 10.7 | 16 |
| Operation at a terrestrial altitude | | max. 2000 | |

Subject to technical changes!

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited