
# IKA

designed for scientists

Dispergierwerkzeuge / Dispersion tools / Outils dispersants

S 10 D - 7G - KS - 65  
S 10 D - 7G - KS - 110

# IKA

designed for scientists

**IKA-Werke GmbH & Co. KG**  
Janke & Kunkel-Straße 10, 79219 Staufen, Germany  
Phone: +49 7633 831-0, Fax: +49 7633 831-98  
eMail: sales@ika.de

**USA**  
IKA Works, Inc.  
Phone: +1 910 452-7059  
eMail: usa@ika.net

**MALAYSIA**  
IKA Works (Asia) Sdn Bhd  
Phone: +60 3 6099-5666  
eMail: sales.lab@ika.my

**JAPAN**  
IKA Japan K.K.  
Phone: +81 6 6730 6781  
eMail: info\_japan@ika.ne.jp

**VIETNAM**  
IKA Vietnam Company Limited  
Phone: +84 28 38202142  
eMail: sales.lab-vietnam@ika.com

**KOREA**  
IKA Korea Ltd.  
Phone: +82 2 2136 6800  
eMail: sales-lab@ika.kr

**CHINA**  
IKA Works Guangzhou  
Phone: +86 20 8222 6771  
eMail: info@ika.cn

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## Safety instructions

- No liquid may get into the drive unit of the disperser.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health and safety regulations accident prevention regulations.
- Do not use any damaged dispersion tools e.g. with hairline cracks or splits.
- Note the operating instructions of the disperser unit.
- Only operate in a vertical working position.

### Product information

The packaging unit encloses 12 or 25 pieces.

The dispersion tool comprises shaft tube and rotor shaft.

Dispersion tools made of plastics are not sterile and not free of DNA or RNA.

The tools are ideal for applications, with which no cross contamination is permissible.

The tools are ideally appropriate for a single - use. However they can be used, depending upon the application, several times too.

In this case please consider the cleaning references.

### Application instruction

**Attention!** Never run the dispersion tool dry. The bearing will be destroyed without cooling by the medium.

Always ensure that the dispersing element does not run dry in some cases should spouts develop.

The fill level of the medium must be between the MIN and MAX marks.

Medium may leak out of the overflow hole around the MAX mark.

The dispersion tool is not suitable for continuous operation.

Dispensing causes the medium to heat up.

Store the dispersing elements in a dry place, in their packaging and at room temperature.

The solid bodies in samples must not be frozen or too hard.

The optimal rotating frequency, usually 20,00-30,000 rpm, must be determined by attempts.

Likewise the optimal dispersion duration must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute.

Longer application times bring no improvements, increase however the sample temperature substantially.

Dispose of the tools after use according to the regulations. Observe the laboratory regulations.

Designation	Application	Max. permissible number of revolution [rpm]	Optimal speed range [rpm]	Working volume [ml]
S10D-7G-KS-65	Liquid/ liquid or soft fabric	30,000	20,000 ... 30,000	1... 20
S10D-7G-KS-110	Liquid/ liquid or soft fabric	30,000	20,000 ... 30,000	1... 40

### Materials and chemical resistance

The shaft tubes are made of polycarbonate (PC), the rotor shafts are made of polysulphone (PSU).

The plastic materials used are approved by FDA (Food and Drug Administration Approval).

The dispersion tool has a good stability against weak acids, chlorides, hypochlorides and many other chemicals.

### Cleaning, disinfection and sterilization

The plastic dispersion tools are not sterile, but you can sterilise them.

	Allowed procedures	Not allowed procedures
Method of sterilisation	Autoclave up to 125 °C with 2 bar <p>Ethylene oxides Formaldehyde (g) Gluteralhyde (2%) Gamma radiation up to 5 MRAD</p>	UV - irradiation Hot-air sterilisation 160–190 °C
Biological disinfection	Alcohol (70 <span> </span> %) Hypochlorides (5 <span> </span> %) Formaldehyde (40 <span> </span> %)	Phenol derivatives Expanded use of guanidine
Radioactive decontamination	Radiacwash Count- off Water/ Ethyl alcohol/ SDS	-

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